

Aptima® Unisex Swab Specimen Collection Kit

Female collection procedure guide


Collection for endocervical swab specimens




Use cleaning swab (white shaft swab with red printing) to remove excess mucus from cervical os and surrounding mucosa. Discard this swab.

Swab specimen collection guide for:

- *Chlamydia trachomatis* (CT)
- *Neisseria gonorrhoeae* (NG)
- *Mycoplasma genitalium* (MG)
- *Trichomonas vaginalis* (TV)



Insert collection swab (blue shaft swab with green printing) into endocervical canal. Gently rotate swab clockwise for 10 to 30 seconds to help ensure adequate sampling. Withdraw swab carefully; avoid any contact with vaginal mucosa.



While holding swab in hand, unscrew the tube cap. Do not spill tube contents. **If the tube contents are spilled, discard and replace with a new Aptima unisex swab transport tube.** Carefully break the swab shaft at the score line against the side of the tube. Discard top portion of swab shaft.



Re-cap swab specimen transport tube tightly.

Specimen Transport and Storage

- After collection, transport and store swab in the unisex specimen transport tube between 2°C to 30°C until tested for Aptima Combo 2® assay, Aptima® Mycoplasma genitalium assay and Aptima® Trichomonas vaginalis assay.
- Specimens must be assayed within 60 days of collection for Aptima Combo 2® assay, Aptima® Mycoplasma genitalium assay and Aptima® Trichomonas vaginalis assay.
- If longer storage is needed, after collection in the Aptima unisex specimen transport tube freeze between -20°C to -70°C for up to 12 months for Aptima Combo 2® assay, up to 90 days for Aptima® Mycoplasma genitalium assay and up to 24 months for the Aptima® Trichomonas vaginalis assay.



Aptima[®] Unisex Swab Specimen Collection Kit


Male collection procedure guide

Collection for male urethral swab specimens

Patient should not have urinated for at least 1 hour prior to specimen collection.


Swab specimen collection guide for:

- *Chlamydia trachomatis* (CT)
- *Neisseria gonorrhoeae* (NG)
- *Mycoplasma genitalium* (MG)



Discard cleaning swab (white shaft with red print on label). The cleaning swab is NOT needed for male specimen collection.

Insert specimen collection swab (blue shaft swab with green printing) 2 cm to 4 cm into urethra. Gently rotate swab clockwise for 2 to 3 seconds in urethra to help ensure adequate sampling. Withdraw swab carefully.



While holding swab in hand, unscrew tube cap. Do not spill tube contents. **If tube contents are spilled, discard and replace with a new Aptima unisex swab transport tube.** Carefully break the swab shaft at the score line against the side of the tube. Discard top portion of swab shaft.



Re-cap swab specimen transport tube tightly.

Specimen Transport and Storage

- After collection, transport and store swab in the unisex specimen transport tube between 2°C to 30°C until tested for Aptima Combo 2[®] assay, Aptima[®] Mycoplasma genitalium assay and Aptima[®] Trichomonas vaginalis assay.
- Specimens must be assayed within 60 days of collection for Aptima Combo 2[®] assay, Aptima[®] Mycoplasma genitalium assay and Aptima[®] Trichomonas vaginalis assay.
- If longer storage is needed, after collection in the Aptima unisex specimen transport tube freeze between -20°C to -70°C for up to 12 months for Aptima Combo 2[®] assay, up to 90 days for Aptima[®] Mycoplasma genitalium assay and up to 24 months for the Aptima[®] Trichomonas vaginalis assay.

Hologic provides this collection procedure guide as a general informational tool only; it is not an affirmative instruction or guarantee of performance. It is the sole responsibility of the clinician/laboratory to read and understand the appropriate package insert and comply with applicable local, state and federal rules and regulations.

hologic.com | diagnostic.solutions@hologic.com | 888.484.4747

DS-01696-001 Rev. 003 © 2019 Hologic, Inc. All rights reserved. Hologic, The Science of Sure, Aptima, Aptima Combo 2 and associated logos are trademarks and/or registered trademarks of Hologic, Inc. and/or its subsidiaries in the United States and/or other countries. This information is intended for medical professionals in the U.S. and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative or write to diagnostic.solutions@hologic.com.

